




LP ITALIANA SPA



CE



Polystyrene spectrophotometer cuvette,
packed in styrofoam boxes of 100 pieces

*Cuvetta in polistirolo per spettrofotometro,
in stativi di polistirolo espanso da 100 pezzi*

Cuve en polystyrène pour spectrophotomètre,
en portoir de polystyrène expansé de 100 pièces

	code	volume ml	dimensions mm	type
1	111117	2-4	10x10x45	macro
2	112117	0.5-2	10x4x45	micro

Polystyrene cup for TECHNICON
type autoanalyzer

*Coppetta in polistirolo per autoanalyzer
tipo TECHNICON*

Godet en polystyrène pour auto-analyseur
type TECHNICON

	code	volume ml	dimensions mm
1	110011	0.5	Ø 14x25
2	110611	2	Ø 14x24
3	110811	4	Ø 17x38



Polystyrene cuvette for SHERWOOD
type coagulometer

*Cuvetta in polistirolo per coagulometro
tipo SHERWOOD*

Godet en polystyrène pour coagulomètre
type SHERWOOD

code	volume ml	dimensions mm
110511	1	Ø 12x27

Polystyrene cup for CENTRIFICHEM
type autoanalyzer

*Coppetta in polistirolo per autoanalyzer
tipo CENTRIFICHEM*

Godet en polystyrène pour auto-analyseur
type CENTRIFICHEM



code	volume ml	dimensions mm
110711	0.25	Ø 14x16

Polystyrene cup for HITACHI
type autoanalyzer

*Coppetta in polistirolo per autoanalyzer
tipo HITACHI*

Godet en polystyrène pour auto-analyseur
type HITACHI



code	volume ml	dimensions mm
110911	3	Ø 17x38


Polypropylene cup for GILFORD 3500
type autoanalyzer

*Coppetta in polipropilene per autoanalyzer
tipo GILFORD 3500*

Godet en polypropylène pour auto-analyseur
type GILFORD 3500



code	volume ml	dimensions mm
110240	0.5	Ø 10x22




1.5 ml polypropylene conical tube, EPPENDORF type
Microprovetta conica da 1,5 ml in polipropilene tipo EPPENDORF
Cupule conique de 1,5 ml en polypropylène type EPPENDORF

	code	graduated graduata graduée	dimensions mm	colour colore couleur
1	111540	■	Ø 10x40	
2	111541	■	Ø 10x40	
3	111543	■	Ø 10x40	
4	111544	■	Ø 10x40	
5	111546	■	Ø 10x40	
6	111549	■	Ø 10x40	
	111542	■	Ø 10x40	Assorted

Conical micro-tube without cap, EPPENDORF type
Microprovetta conica senza tappo, tipo EPPENDORF
Cupule conique sans bouchon, type EPPENDORF

	code	material	volume ml	dimensions mm
1	111610	PS	1.5	Ø 10x40
2	111640	PP	1.5	Ø 10x40



Polyethylene conical micro-tube with cap, BECKMAN type

Microprovetta conica in polietilene con tappo, tipo BECKMAN

Cupule conique en polyéthylène avec bouchon, type BECKMAN



code	volume ml	dimensions mm
110540	0.5	Ø 6x44

HDPE conical micro-tube with cap, COBAS type
Microprovetta conica in HDPE con tappo, tipo COBAS
 Cupule conique en HDPE avec bouchon, type COBAS

	code	volume ml	dimensions mm	colour colore couleur
1	110340	0.75	Ø 8x35	
2	110343	0.75	Ø 8x35	
3	110344	0.75	Ø 8x35	
4	110345	0.75	Ø 8x35	
5	110346	0.75	Ø 8x35	



Polypropylene conical micro-tube with cap,
 Vitatron Akes type

*Microprovetta conica in polipropilene con tappo,
 tipo Vitatron Akes*

Cupule conique en polypropylène avec bouchon,
 type Vitatron Akes

code	volume ml	dimensions mm
110549	0.5	Ø 8x31



Cuvette segment for COBAS MIRA type
 autoanalyzer
*Cuvetta segmento per autoanalyzer tipo
 COBAS MIRA*
 Segment de cuves pour auto-analyseur type
 COBAS MIRA

	code	
1	110357	
2	111357	15 pcs. rack

