



LP ITALIANA SPA



Anticoagulant tubes
Provette con anticoagulante
Tubes avec anticoagulant

Sodium iodoacetate



KF + Na₂EDTA



K₂EDTA



K₃EDTA



In 114040, polypropylene, Ø 12x51 mm



113307

for 2.5 ml
of blood
per 2,5 ml
di sangue
pour 2,5 ml
de sang

113317

for 2.5 ml
of blood
per 2,5 ml
di sangue
pour 2,5 ml
de sang

113337

for 2.5 ml
of blood
per 2,5 ml
di sangue
pour 2,5 ml
de sang

113397

for 2.5 ml
of blood
per 2,5 ml
di sangue
pour 2,5 ml
de sang

In 110440, polypropylene, Ø 12x75 mm



110407

for 2.5 ml
of blood
per 2,5 ml
di sangue
pour 2,5 ml
de sang

110417

for 2.5 ml
of blood
per 2,5 ml
di sangue
pour 2,5 ml
de sang

110437

for 2.5 ml
of blood
per 2,5 ml
di sangue
pour 2,5 ml
de sang

110487

for 2.5 ml
of blood
per 2,5 ml
di sangue
pour 2,5 ml
de sang

In 115040, polypropylene, Ø 16x60 mm



114307

for 2.5 ml
of blood
per 2,5 ml
di sangue
pour 2,5 ml
de sang

114317

for 2.5 ml
of blood
per 2,5 ml
di sangue
pour 2,5 ml
de sang


114337

for 2.5 ml
of blood
per 2,5 ml
di sangue
pour 2,5 ml
de sang

114397

for 2.5 ml
of blood
per 2,5 ml
di sangue
pour 2,5 ml
de sang

In 113040, polypropylene, Ø 12x92 mm



115307

for 2.5 ml
of blood
per 2,5 ml
di sangue
pour 2,5 ml
de sang

115317

for 2.5 ml
of blood
per 2,5 ml
di sangue
pour 2,5 ml
de sang

115337

for 2.5 ml
of blood
per 2,5 ml
di sangue
pour 2,5 ml
de sang

115397

for 2.5 ml
of blood
per 2,5 ml
di sangue
pour 2,5 ml
de sang

In 110440, polypropylene, Ø 12x75 mm



111407

for 5 ml
of blood
per 5 ml
di sangue
pour 5 ml
de sang

111417

for 5 ml
of blood
per 5 ml
di sangue
pour 5 ml
de sang

111437

for 5 ml
of blood
per 5 ml
di sangue
pour 5 ml
de sang

111497

for 5 ml
of blood
per 5 ml
di sangue
pour 5 ml
de sang

In 115040, polypropylene, Ø 16x60 mm



116307

for 5 ml
of blood
per 5 ml
di sangue
pour 5 ml
de sang

116317

for 5 ml
of blood
per 5 ml
di sangue
pour 5 ml
de sang


116337

for 5 ml
of blood
per 5 ml
di sangue
pour 5 ml
de sang

116397

for 5 ml
of blood
per 5 ml
di sangue
pour 5 ml
de sang

In 113040, polypropylene, Ø 12x92 mm



117307

for 5 ml
of blood
per 5 ml
di sangue
pour 5 ml
de sang

117317

for 5 ml
of blood
per 5 ml
di sangue
pour 5 ml
de sang


117337

for 5 ml
of blood
per 5 ml
di sangue
pour 5 ml
de sang

117397

for 5 ml
of blood
per 5 ml
di sangue
pour 5 ml
de sang


In 119040, polypropylene, Ø 16x75 mm



117497

Hemalog, for
5 ml of blood
per 5 ml
di sangue
pour 5 ml
de sang

In 111040, polypropylene, Ø 16x100 mm



118307

for 5 ml
of blood
per 5 ml
di sangue
pour 5 ml
de sang

118317

for 5 ml
of blood
per 5 ml
di sangue
pour 5 ml
de sang

118337

for 5 ml
of blood
per 5 ml
di sangue
pour 5 ml
de sang

118397

for 5 ml
of blood
per 5 ml
di sangue
pour 5 ml
de sang

Sodium Heparin	Lithium Heparin	Sodium Citrate	Sodium Citrate	Sodium Citrate	Sodium Citrate	ACD formule A	ACD formule A
							
113327 for 2.5 ml of blood <i>per 2,5 ml di sangue</i> pour 2,5 ml de sang	113357 for 2.5 ml of blood <i>per 2,5 ml di sangue</i> pour 2,5 ml de sang		113577 0.25 ml	113877 0.4 ml	113477 0.5 ml		
110427 for 2.5 ml of blood <i>per 2,5 ml di sangue</i> pour 2,5 ml de sang	110457 for 2.5 ml of blood <i>per 2,5 ml di sangue</i> pour 2,5 ml de sang			110477 0.4 ml			
114327 for 2.5 ml of blood <i>per 2,5 ml di sangue</i> pour 2,5 ml de sang	114357 for 2.5 ml of blood <i>per 2,5 ml di sangue</i> pour 2,5 ml de sang	114867 0.4 ml	114577 0.25 ml	114877 0.4 ml	114477 0.5 ml		
115327 for 2.5 ml of blood <i>per 2,5 ml di sangue</i> pour 2,5 ml de sang	115357 for 2.5 ml of blood <i>per 2,5 ml di sangue</i> pour 2,5 ml de sang	115867 0.4 ml	115577 0.25 ml	115877 0.4 ml	115477 0.5 ml		
111427 for 5 ml of blood <i>per 5 ml di sangue</i> pour 5 ml de sang	111457 for 5 ml of blood <i>per 5 ml di sangue</i> pour 5 ml de sang				111477 0.5 ml		
116327 for 5 ml of blood <i>per 5 ml di sangue</i> pour 5 ml de sang	116357 for 5 ml of blood <i>per 5 ml di sangue</i> pour 5 ml de sang	114367 0.5 ml			114377 1 ml		
117327 for 5 ml of blood <i>per 5 ml di sangue</i> pour 5 ml de sang	117357 for 5 ml of blood <i>per 5 ml di sangue</i> pour 5 ml de sang	115367 0.5 ml			115377 1 ml	117387 0.75 ml	117287 1.5 ml
118327 for 10 ml of blood <i>per 10 ml di sangue</i> pour 10 ml de sang	118357 for 10 ml of blood <i>per 10 ml di sangue</i> pour 10 ml de sang	118367 0.5 ml			118377 1 ml	118387 0.75 ml	118287 1.5 ml



Autoplast® K₃EDTA anticoagulant for 3 ml of blood with pierceable stopper and label
Autoplast® anticoagulante K₃EDTA per 3 ml di sangue con tappo perforabile ed etichetta
Autoplast® K₃EDTA pour 3 ml de sang avec bouchon perforable et étiquette

code	material	volume ml	dimensions mm
110497	PP	5	Ø 12x75

0.35 ml sodium citrate for ACL 2000 type analyzer for 3.15 ml of blood
0,35 sodio citrato per autoanalyzer tipo ACL 2000 per 3,15 ml di sangue
0,35 ml citrate de sodium pour analyseur type ACL 2000 pour 3,15 ml de sang



code	material	volume ml	dimensions mm
110467	PP	5	Ø 12x75




Miniplast® conical mini test tube, Ø 12x51 mm, complete with sterilized capillary drawing system
Miniplast® miniprovetta conica, Ø 12x51 mm, completa di unità di prelievo sterile
Miniplast® mini-tube, Ø 12x51 mm, avec système de prélèvement stérile

	code	anticoagulant	blood ml
1	241010	no	1.5
2	242097	K ₃ EDTA	0.6
3	242057	Lithium Heparin	0.6




Separgel® Elleplast test tube with separating silicone gel and blood clot accelerator
Separgel® provetta in Elleplast con gel separatore e acceleratore della coagulazione
Separgel® tube en Elleplast avec gel séparateur et accélérateur de la coagulation

	code	volume ml	dimensions mm	label etichetta étiquette	blood clot accelerator acceleratore coagulazione accélérateur coagulation
1	221050	10	Ø 16x100		
2	221150	10	Ø 16x100	■	
3	228050	5	Ø 12x100		
4	228150	5	Ø 12x100	■	
5	223050	5	Ø 12x88		
6	223150	5	Ø 12x88	■	
7	231059	10	Ø 16x100		■
8	231159	10	Ø 16x100	■	■
9	238059	5	Ø 12x100		■
10	238159	5	Ø 12x100	■	■
11	233059	5	Ø 12x88		■
12	233159	5	Ø 12x88	■	■



Separplast® test tube with polystyrene separating granules
Separplast® provetta con granuli separatori di polistirolo
Separplast® tube avec granules de séparation en polystyrène

	code	material	volume ml	dimensions mm	label etichetta étiquette
1	211042	PP	10	Ø 16x100	
2	211142	PP	10	Ø 16x100	■
3	211012	PS	10	Ø 16x100	
4	211112	PS	10	Ø 16x100	■
5	211052	Elleplast	10	Ø 16x100	
6	211152	Elleplast	10	Ø 16x100	■
7	218052	Elleplast	5	Ø 12x100	
8	218152	Elleplast	5	Ø 12x100	■
9	213052	Elleplast	5	Ø 12x88	
10	213152	Elleplast	5	Ø 12x88	■



Separplast plus® test tube with polystyrene separating granules covered with clot accelerator

Separplast plus® provetta con granuli separatori di polistirolo ricoperti da un acceleratore della coagulazione

Separplast plus® tube avec granules de séparation en polystyrène couverts d'un accélérateur de la coagulation

	code	material	volume ml	dimensions mm	label etichetta étiquette
1	251045	PP	10	Ø 16x100	
2	251145	PP	10	Ø 16x100	■
3	251015	PS	10	Ø 16x100	
4	251115	PS	10	Ø 16x100	■
5	251055	Elleplast	10	Ø 16x100	
6	251155	Elleplast	10	Ø 16x100	■
7	258055	Elleplast	5	Ø 12x100	
8	258155	Elleplast	5	Ø 12x100	■
9	253055	Elleplast	5	Ø 12x88	
10	253155	Elleplast	5	Ø 12x88	■