



LP ITALIANA SPA





High grip and closed bottom
polyethylene stopper

*Tappo in polietilene a presa alta e
chiuso sul fondo*


Bouchon en polyéthylène à bord haut
et fermé au fond


code	dimensions mm	colour colore couleur
116130	Ø 16	
116131	Ø 16	
116132	Ø 16	
116133	Ø 16	
116134	Ø 16	
116135	Ø 16	
116136	Ø 16	
116141	Ø 16	
116142	Ø 16	
116143	Ø 16	
116145	Ø 16	
116146	Ø 16	
116148	Ø 16	

High grip and closed bottom
polyethylene stopper


*Tappo in polietilene a presa alta e
chiuso sul fondo*

Bouchon en polyéthylène à bord haut
et fermé au fond

code	dimensions mm	colour colore couleur
118130	Ø 12	
118131	Ø 12	
118132	Ø 12	
118133	Ø 12	
118134	Ø 12	
118135	Ø 12	
118136	Ø 12	
118141	Ø 12	
118142	Ø 12	
118143	Ø 12	
118145	Ø 12	
118146	Ø 12	
118148	Ø 12	



Polyethylene winged stopper
Tappo alettato in polietilene
Bouchon à ailettes en polyéthylène

code	dimensions mm	colour colore couleur
117130	Ø 16	
117131	Ø 16	
117132	Ø 16	
117133	Ø 16	
117134	Ø 16	
117136	Ø 16	



Polyethylene winged stopper
Tappo alettato in polietilene
Bouchon à ailettes en polyéthylène

code	dimensions mm	colour colore couleur
119130	Ø 11-13	
119131	Ø 11-13	
119132	Ø 11-13	
119133	Ø 11-13	
119134	Ø 11-13	
119136	Ø 11-13	





Red rubber conical stopper
Tappo conico in gomma rossa
Bouchon conique en caoutchouc rouge

code	upper mm superiore mm supérieur mm	lower mm inferiore mm inférieur mm	height mm altezza mm hauteur mm
291200	Ø 9	Ø 6	18
291201	Ø 12	Ø 9	18
291202	Ø 13	Ø 10	18
291203	Ø 14	Ø 11	18
291204	Ø 16	Ø 12	20
291205	Ø 18	Ø 13	22
291206	Ø 20	Ø 14	24
291207	Ø 23	Ø 16	26
291208	Ø 26	Ø 19	28
291209	Ø 30	Ø 21	30
291210	Ø 32	Ø 25	34
291211	Ø 37	Ø 28	38
291212	Ø 42	Ø 32	42
291213	Ø 45	Ø 34	45
291214	Ø 48	Ø 37	50
291215	Ø 54	Ø 42	53
291216	Ø 60	Ø 43	60
291217	Ø 65	Ø 53	70

Silicone conical stopper
Tappo conico in silicone
Bouchon conique en silicone

code	upper mm superiore mm supérieur mm	lower mm inferiore mm inférieur mm	height mm altezza mm hauteur mm
291220	Ø 9	Ø 6	18
291221	Ø 12	Ø 9	18
291222	Ø 13	Ø 10	18
291223	Ø 14	Ø 11	18
291224	Ø 16	Ø 12	20
291225	Ø 18	Ø 13	22
291226	Ø 20	Ø 14	24
291227	Ø 23	Ø 16	26
291228	Ø 26	Ø 19	28
291229	Ø 30	Ø 21	30
291230	Ø 32	Ø 25	34
291231	Ø 37	Ø 28	38
291232	Ø 42	Ø 32	42
291233	Ø 45	Ø 34	45
291234	Ø 48	Ø 37	50
291235	Ø 54	Ø 42	53
291236	Ø 60	Ø 43	60



Polyethylene stopper for test tubes
Tappo in polietilene per provette
Bouchon en polyéthylène pour tubes à essai

code	internal mm interno mm intérieur mm	test tube mm provetta mm tube mm
290350	Ø 10.5	Ø 12
290351	Ø 14	Ø 16
290353	Ø 16	Ø 18
290354	Ø 18	Ø 20
290355	Ø 22	Ø 24
290356	Ø 28	Ø 30
290357	Ø 33	Ø 35
290358	Ø 38	Ø 40
290359	Ø 43	Ø 45
290360	Ø 48	Ø 50



Polyethylene normalized tapered stopper
Tappo conico normalizzato, in polietilene
Bouchon conique normalisé, en polyéthylène



code	NS
291080	7/16
291081	10/19
291082	12/21
291083	14/23
291084	19/26
291085	24/29
291086	29/32
291087	34/35
291088	45/40
291089	60/50